

CELLFREE[®] Human sIL-2R ELISA Kit

EH2IL2R

1352.5

Number	Description
EH2IL2R	<p>CELLFREE[®] Human Soluble Interleukin-2 Receptor (sIL-2R) ELISA, sufficient reagents for 96 determinations</p> <p>Kit contents:</p> <p>Anti-Human sIL-2R Pre-coated 96-well Strip Plate, 1 each</p> <p>Lyophilized Prediluted Recombinant Human sIL-2R Standards, 2 vials</p> <p>Lyophilized Recombinant Human sIL-2R Controls, 2 vials</p> <p>Sample Diluent, 1 vial</p> <p>Conjugate Reagent, 1 vial</p> <p>20X Wash Buffer, 1 vial</p> <p>TMB Substrate, 1 vial</p> <p>Stop Solution, 13 ml, 0.16 M sulfuric acid</p> <p>Adhesive Plate Covers, 2 each</p>

For research use only – not for use in diagnostic procedures.

Storage: For maximum stability, store in a non-defrosting -20°C freezer and refer to the expiration date for frozen storage on the label. Alternatively, store at 2-8°C and refer to the expiration date for refrigerated storage. Once thawed, store at 4°C until the expiration date for refrigerated storage. Kit is shipped on dry ice.









Introduction

The Thermo Scientific CELLFREE Human Soluble Interleukin-2 Receptor (sIL-2R) ELISA is an *in vitro* enzyme-linked immunosorbent assay for the quantitative measurement of human sIL-2R in serum; EDTA, sodium citrate and heparin plasma and culture supernatants.

Table of Contents

Introduction	1
Procedure Summary.....	2
Additional Materials Required.....	2
Precautions.....	2
Sample Preparation	3
Reagent Preparation.....	3
Assay Procedure	4
Performance Characteristics	5
References	6
Data Templates	7

Procedure Summary

- | | | | |
|---|---|---|--|
| 
1. Add 50 µl of Standards or samples to each well in duplicate. | 
2. Add 50 µl of Conjugate Reagent to each well. | 
3. Cover plate and incubate at room temperature (20-25°C) for 2 hours on a rotator set at 150 (±10) RPM. | 
4. Wash plate THREE times. |
| 
5. Add 100 µl of TMB Substrate to each well. | 
6. Incubate plate in the dark at room temperature for 30 minutes. | 
7. Stop reaction by adding 100 µl of Stop Solution to each well. | 
8. Measure absorbance on a plate reader at 450 nm or 450 minus 550 nm and calculate results. |

Additional Materials Required

- Precision pipettors with disposable plastic tips to deliver 5-1,000 µl
- Plastic pipettes to deliver 5-15 ml
- Ultrapure water
- A glass or plastic one-liter container to prepare Wash Buffer
- A squirt wash bottle or an automated 96-well plate washer
- 1.5 ml polypropylene or polyethylene tubes to prepare standards – do not use polystyrene, polycarbonate or glass tubes
- Disposable reagent reservoirs, 3 each
- A rotator set at 150 (±10) RPM
- A standard ELISA reader for measuring absorbance at 450 nm and 550 nm. If a 550 nm filter is not available, the absorbance can be measured at 450 nm only. Refer to the instruction manual supplied with the instrument being used.
- Graph paper or a computerized curve-fitting statistical software package

Precautions

- **All samples and reagents must be at room temperature (20-25°C) before use in the assay.**
- Review all instructions carefully and verify all components against the Kit Contents list (page 1) before beginning.
- Sodium azide inactivates the Conjugate Reagent. **DO NOT USE SAMPLES CONTAINING SODIUM AZIDE IN THIS ASSAY.**
- Do not use a water bath to thaw samples. Thaw samples at room temperature.
- If using a multichannel pipettor, always use a new disposable reagent reservoir.
- Use new disposable pipette tips for each transfer to avoid cross-contamination.
- Use a new adhesive plate cover for each incubation step.
- Once reagents have been added to the plate, take care **NOT** to let plate **DRY** at any time during the assay.
- Avoid microbial contamination of reagents.
- Vigorous plate washing is essential.
- Avoid exposing reagents to excessive heat or light during storage and incubation.
- Do not mix reagents from different kit lots. Discard unused ELISA components after assay completion.

- Do not use glass pipettes to measure TMB Substrate. Take care not to contaminate the solution. If the solution is blue before use, DO NOT USE IT.
- Crystals may form in the 20X Wash Buffer during storage. Dissolve crystals before diluting by warming in a 36-38°C water bath and mixing.
- Individual components may contain antibiotics and preservatives. Wear gloves while performing the assay to avoid contact with samples and reagents. Please follow proper disposal procedures.

Sample Preparation

Sample Handling

- Serum; EDTA, sodium citrate and heparin plasma; or culture supernatants may be tested in this ELISA.
- 50 µl per well of serum, plasma or culture supernatant are required.
- Store samples to be assayed within 24 hours at 2-8°C. For long-term storage, aliquot and freeze samples at -70°C.
- Avoid repeated freeze-thaw cycles when storing samples. Mix samples by gently inverting tubes.
- Test samples and standards must be assayed in duplicate each time the assay is performed.
- Gradually equilibrate samples to room temperature before beginning assay. Do not use a heated water bath to thaw or warm samples.
- If samples are clotted, grossly hemolyzed, lipemic or microbially contaminated, or if there is any question about the integrity of a sample, make a note on the template and interpret results with caution.

Sample Dilution

- If the human sIL-2R concentration possibly exceeds the highest point of the standard curve (i.e., ~6,500 U/ml), prepare one or more 10-fold dilutions of the test sample. When testing culture supernatants, prepare the serial dilutions using your culture medium. When testing serum or plasma prepare the serial dilutions using the Sample Diluent provided. For example, a 10-fold dilution is prepared by adding 50 µl of test sample to 450 µl of appropriate diluent. Mix thoroughly between dilutions before assaying.

Reagent Preparation

Look for (PP) throughout these instructions for procedural differences when using partial plates.

Wash Buffer

1. Label a clean glass or plastic one-liter container "Wash Buffer."
2. Add the entire contents of the 20X Wash Buffer bottle to the one-liter container and dilute to a final volume of 1.0 liter with ultrapure water. Mix thoroughly.

(PP) When using partial plates, store the reconstituted Wash Buffer at 2-8°C.

Note: Crystals may form in the 20X Wash Buffer during storage. Dissolve crystals before dilution by warming in a 36-38°C water bath and mixing.

Note: Wash buffer must be at room temperature before use in the assay. Do not use Wash Buffer if it becomes visibly contaminated during storage.


Standards and Controls

1. Reconstitute one vial of Lyophilized Standard per partial plate. See the vial label for the reconstitution volume.
2. Reconstitute Lyophilized Controls with 0.5 ml of ultrapure water.
3. Let stand for 10 minutes to dissolve completely. Gently swirl vial to mix.

Note: Reconstituted standards and controls may be stored at 2-8°C for up to 7 days. If reconstituted standards are not used within 7 days, store at -70°C until kit expiration date.

Note: Unit values for the controls vary slightly from lot to lot. See concentrations stated on each vial label to determine control levels.

4. Label six tubes, one for each standard curve point: 6,500 U/ml, 3,250 U/ml, 1,625 U/ml, 812.5 U/ml, 406.3 U/ml and 0 U/ml, then prepare 1:2 serial dilutions for the standard curve as follows:
5. Pipette 200 µl of appropriate diluent into each tube.
6. Pipette 200 µl of the reconstituted standard into the first tube (i.e., 6,500 U/ml) and mix.
7. Pipette 200 µl of this dilution into the second tube (i.e., 3,250 U/ml) and mix.
8. Repeat the serial dilutions (using 200 µl) three more times to complete the standards.



Assay Procedure

A. Standard, Sample and Conjugate Reagent Incubation

- **(PP)** Determine number of strips required. Leave these strips in the plate frame. Tightly seal unused strips in the provided foil pouch with desiccant and store at 2-8°C. After completing assay, retain plate frame for second partial plate. When using the second partial plate, place strips securely in the plate frame.
 - Use the Data Template provided to record the locations of the zero standard (blank or negative control), human sIL-2R standards, controls and test samples. Perform five standard points and one blank in duplicate with each series of unknown samples. Two controls are provided with this assay.
 - If using a multichannel pipettor, use a new reagent reservoir to add the Conjugate Reagent. Remove from the vial only the amount required for the number of strips being used. Take care not to touch the samples in wells with the pipette tip when adding the Conjugate Reagent.
 - **(PP)** Do not combine leftover Conjugate Reagent with that reserved for the other partial plate.
1. Add 50 µl of reconstituted standards, controls or test samples to each well in duplicate. Mix well by gently tapping the plate several times.

Note: If the human sIL-2R concentration in any test sample possibly exceeds the highest point on the standard curve, 6,500 U/ml, see Sample Preparation – Sample Dilution section.

2. Add 50 µl of Conjugate Reagent to each well.
3. Carefully cover plate with an adhesive plate cover. Ensure all edges and strips are tightly sealed by running your thumb over edges and down each strip. Incubate for two (2) hours at room temperature, 20-25°C, on a rotator set at 150 (±10) RPM.
4. Carefully remove adhesive plate cover. Wash plate THREE times with Wash Buffer as described in the Plate Washing section (section B).

B. Plate Washing

1. Gently squeeze the long sides of plate frame before washing to ensure all strips securely remain in the frame.
2. Discard plate contents. Use a squirt wash bottle to vigorously fill each well completely with Wash Buffer, then discard plate contents. Repeat procedure two additional times for a total of THREE washes. Blot plate onto paper towels or other absorbent material.

Note: For automated washing, aspirate all wells and wash THREE times with Wash Buffer, overfilling wells with Wash Buffer. Blot plate onto paper towels or other absorbent material.

C. Substrate Incubation and Stop Step

- Use new disposable reagent reservoirs when adding TMB Substrate Solution and Stop Solution.
 - Dispense from bottle ONLY amount required, 100 µl per well, for the number of wells being used. Do not use a glass pipette to measure the TMB Substrate Solution.
 - **(PP)** Do not combine leftover substrate with that reserved for the second partial plate. Take care not to contaminate remaining TMB Substrate Solution.
1. Pipette 100 µl of TMB Substrate Solution into each well.
 2. Allow enzymatic color reaction to develop at room temperature in the dark for 30 minutes. Do not cover plate with aluminum foil or a plate sealer. The substrate reaction yields a blue solution that turns yellow when Stop Solution is added.
 3. After 30 minutes, stop the reaction by adding 100 µl of Stop Solution to each well.

D. Absorbance Measurement


- **Evaluate the plate within 30 minutes of stopping the reaction.** Measure absorbance on an ELISA plate reader set at 450 nm and 550 nm. Subtract 550 nm values from 450 nm values to correct for optical imperfections in the microplate. If an absorbance at 550 nm is unattainable, measure at 450 nm only.

Note: When the 550 nm measurement is omitted, absorbance values will be higher.

E. Calculation of Results

- Generate the standard curve by plotting the average absorbance obtained for each Standard concentration on the vertical (Y) axis vs. the corresponding human sIL-2R concentration (pg/ml) on the horizontal (X) axis.
- Calculate results using graph paper or curve-fitting statistical software. The human sIL-2R amount in each sample is determined by interpolating from the absorbance value (Y axis) to human sIL-2R concentration (X axis) using the standard curve.
- If the test sample was diluted, multiply the interpolated value obtained from the standard curve by the dilution factor to calculate U/ml of human sIL-2R in the sample.
- Absorbance values obtained for duplicates should be within 10% of the mean value. Carefully consider duplicate values that differ from the mean by greater than 10%.

Example Standard Curve



Performance Characteristics

Sensitivity: < 24 U/ml

The sensitivity or lower limit of detection (LLD)¹ was determined by assaying replicates of zero and the standard curve. The mean signal of zero + 2 standard deviations read in dose from the standard curve is the LLD. This value is the smallest dose that is not zero with 95% confidence.

Specificity:

This ELISA is specific for the measurement of natural and recombinant human sIL-2R. It does not cross react with human IL-1α, IL-1β, IL-2, IL-3, IL-4, IL-5, IL-6, IL-7, IL-8, IL-9, IL-10, IL-11, IL-12, IL-13, IL-15, GM-CSF, IFNγ, IFNα, TNFβ, or mouse IL-2.

Calibration: The standards in this ELISA are calibrated to an internal sIL-2R reference preparation of supernatant from a transfected cell line that secretes recombinant soluble IL-2R. We define 1,000 units as the amount of sIL-2R present in 1 ml of an internal reference preparation. Units/ml can be converted to sIL-2R pg/ml using the following formula:

$$pg/ml = \frac{value(units/ml)}{0.113}$$

Expected Values: Serum and plasma samples were collected from healthy individuals and evaluated in this assay. The levels of sIL-2R found in these samples are indicated in Table 1.

Table 1. Human sIL-2R levels from apparently healthy individuals.

<u>Sample</u>	<u>Average</u>	<u>Range</u>
Serum (n=18)	521 U/ml	269-1,116 U/ml
Plasma (n=20)	642 U/ml	342-1,521 U/ml

Assay Range: 0 to 6,500 U/ml

Unit values of sIL-2R in Prediluted Standards will vary slightly depending on the kit lot.

Reproducibility:



Intra-assay CV: < 10%

Inter-assay CV: < 10%

Precision: Typical intra-assay coefficient of variation is plotted against sIL-2R concentration (U/ml). The points represent samples and controls evaluated in duplicate.

Recovery: Recovery was determined by spiking 1,56487 U/ml natural human sIL-2R into 19 plasma samples and 20 serum samples collected from apparently healthy individuals, and a Sample Diluent control buffer. Endogenous sIL-2R levels were subtracted from the spiked values. Mean serum recovery was 121% and mean plasma recovery was 126%.

Linearity of Dilution: Dilution linearity was determined by serially diluting nine different positive samples. The dilutions were evaluated in the ELISA and the “observed” doses are plotted against the “expected” doses.



Reference

1. *Immunoassay: A Practical Guide*, Chan and Perlstein, Eds., 1987, Academic Press: New York, p71.

Data Templates

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current versions of product instructions are available at www.thermo.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2009 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.